20-125-9496

May 04, 2020
RECEIVED DATE
May 02, 2020

SEND TO **35432**



PAGE 1/2

| SSUE DATE | May 04, 2020

13611 B Street • Omaha, Nebraska 68144-3693 • (402) 334-7770 www.midwestlabs.com

REPORT OF ANALYSIS

For: (35432) OPEN FARM
Cat-OF Whitefish & Green Lentils Kibble
12120FCF

	Level Found		Reporting			Analyst-	Verified-
Analysis	As Receiv	ved Uni	its Lim	nit	Method	Date	Date
Sample ID: 12120FCF-Micros	Lab Number: 13429228	Date Sample	d: 2020-04-	30			
Listeria	negati	ive org/	125g	1	RapidChek/AOAC RI 020401	clh4-2020/05/04	jzh4-2020/05/04
E. coli (generic)	n	.d. cf	u/g	10	AOAC 991.14	jzh4-2020/05/04	arj0-2020/05/04
Salmonella	negati	ive org/	′375g	1	RapidChek/AOAC RI 030301; AFNOR SDI 34/01-04/10	clh4-2020/05/04	jzh4-2020/05/04

All results are reported on an AS RECEIVED basis., n.d. = not detected, cfu = colony forming unit

For questions please contact:

Sara Sterkel-Colombo Account Manager

scolombo@midwestlabs.com (402)590-2983

REPORT NUMBER

20-125-9496

May 04, 2020 RECEIVED DATE May 02, 2020 35432



PAGE 2/2

May 04, 2020

13611 B Street • Omaha, Nebraska 68144-3693 • (402) 334-7770 www.midwestlabs.com

REPORT OF ANALYSIS

For: (35432) OPEN FARM Cat-OF Whitefish & Green Lentils Kibble 12120FCF

Detailed Method Description(s)

E. coli and Total Coliform using 3M Petrifilm

Sample analysis follows MWL MI 292 which is based on AOAC 2018.13. A representative sample is obtained and added to phosphate buffer. Aliquots of the sample are withdrawn and placed on Petrifilm plates. The plates are incubated for 18 to 24 hours. After incubation, the plates are counted to determine the number of generic E. coli and total coliforms present. The color of the colony and the presence of gas differentiate a generic coliform from E. coli. The levels are reported as colony forming units (cfu).

Salmonella - Lateral Flow

Samples are analyzed following MWL MI 195 which is based on the RapidChek Select Salmonella User Guide. A representative sample is obtained and combined with a selective media and allowed to incubate. After incubation, a test strip is used for Salmonella determination. Results are reported as negative or presumptive positive.

Listeria Lateral Flow

Samples are analyzed following MWL MI 194 which is based on the RapidChek Listeria User Guide. A representative sample is obtained and combined with a selective growth media. It is incubated for 40-48 hours. After incubation, an aliquot is heated for 10 minutes, and a test strip for Listeria detection is used. Results are reported as negative or presumptive positive. This procedure does not speciate Listeria.